

Express Installation Instructions



Note: These instructions are provided for the experienced Micromedex® Healthcare Series customer. For detailed, step-by-step installation instructions, reference or print the **Install.pdf** file located in the \document folder on the CD-ROM labeled, "Install".

The Healthcare Series data and software is contained on five CD-ROMs; four data CDs (Disc A, B, C, and D) and one install CD. Each quarterly update is a complete install; meaning that an update installation can and should be installed into the same directory and folders as the previous release.

Standalone Workstation Installation

1. Insert the CD labeled "Install" into the local CD-ROM drive.
2. Run the **setup.exe** located on the Install disc.
3. Accept the default options and parameters whenever possible to minimize update issues. The default data path for disc A is: `C:\mdxw\hcs32\hcsa` (...hcsb for disc B, ...hcsd disc C, and ...hcsd disc D).
4. **If** you are prompted for the Customer ID and Activation Code, insert the diskette sent with this volume, (see Note below) or type the codes which can be found on the software packing slip.



Note: The floppy diskette is volume-specific; meaning that a floppy sent with Volume 100 can only be used to install the Volume 100 data and software.

5. When asked which data discs you wish to copy, select all four if drive space is available.
6. If you chose to copy the data discs in step 5, you are prompted to insert the first data disc. The next screen asks for the drive letter of the CD-ROM drive where the data CD has been inserted. For example: Type, `D:\` if that is the drive letter where the Healthcare Series data CD has been inserted. Step 6 is repeated for *each* data CD you chose to copy in step 5.

Network Installation

The network installation should be performed *from* a workstation *to* the server, using the same path that will be used by your Healthcare Series users. Both UNC and mapped drives are supported:

- UNC Example: `\\<netpath>\apps\mdxw\hcs32`
- Mapped Drive Example: `<drive-letter>:\apps\mdxw\hcs32`

1. Insert the CD labeled "Install" into the local CD-ROM drive.
2. Run the **setup.exe** located on the Install disc.
3. Accept the default options and parameters whenever possible to minimize update issues. The default data path for disc A is: `<netpath>\mdxw\hcs32\hcsa` (...hcsb for disc B, ...hcsd, and ...hcsd disc D).
4. Customer ID & Activation Code: Same as Standalone Step 4 instructions.
5. Copy the data CDs: Same as Step 5 and 6 of Standalone instructions, however, if data transfer rates from the workstation to the server are a concern, you may want to run the copy process manually. To copy manually, do not select any of the data CDs in this step. When ready to copy manually, go to the server and copy the data CDs to the data paths specified in step 3.
Run the client/workstation installation on any *new* workstations after running the network update installation. *Existing* workstations normally do not require installation after an update installation.

Network Client/Workstation Installation

1. Run: `<netpath>\mdxw\hcs32\setup\setup.exe /s` from a command line at the workstation to perform a "silent" (accepts all defaults) installation.